

PRODUCT INFORMATION

- 1 Flat lateral design
- 2 Smooth radius transition from proximal to distal portion of the nail
- 3 Flutes for stress reduced insertion
- 4 Distal diameter: 9, 10, 11 and 12 mm
- 5 CCD angle: 125° and 130°
- 6 Nail length: 170 mm, 200 mm, 240 mm left & right, 300 – 480 mm left & right



HERACLES PROXIMAL FEMUR NAIL



7S MEDICAL AG

Länggasse 4
6208 Oberkirch

+41 41 925 66 60
+41 41 588 01 99

info@7smedical.com
www.7smedical.com

7Smedical 
Swiss Experts in Orthopedics.

HERACLES PROXIMAL FEMUR NAIL



One System Providing Patient Specific Treatment Options

- Blade with three anchors for stable fixation in osteoporotic bone.
- Lag Screw with proven thread design.

Straightforward and Accurate Insertion


One piece radiolucent aiming arm:

- Less weight.
- Increased width to compensate excessive soft tissue.
- Aiming support possibilities.




HERACLES Long Nail

Smooth Nail Insertion to Preserve Reduction




Profile conventional nail.




Flat lateral design.

Just Lock with TwinLock



- Preassembled locking mechanism.
- Controlled and limited gliding.
- Locking of rotation.



Unique End Cap for the option to stop gliding for certain indications.